




[www.acmastechnologies.com](http://www.acmastechnologies.com) | [www.acmas.in](http://www.acmas.in)

# LAB WILLEY GRINDER (ARTHUR H. THOMAS TYPE)



**ACMAS Technologies Inc.**

CE | An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA  
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603  
Email: [acmastechnologies@gmail.com](mailto:acmastechnologies@gmail.com) | [obromax.delhi@gmail.com](mailto:obromax.delhi@gmail.com)  
[www.acmastechnologies.com](http://www.acmastechnologies.com) | [www.acmas.in](http://www.acmas.in)

**ACMAS Technologies Inc.**


CE | An ISO 9001:2000 Company



# LAB WILLEY GRINDER (ARTHUR H. THOMAS TYPE)

## Optional Accessories

a) Automatic 0-60 minute " German" timer



# LAB WILLEY GRINDER (ARTHUR H. THOMAS TYPE)

The chamber is fitted with a hopper provided with a cover. Three delivery S.S. sieves, of 20, 40, 60 mesh and plastic receiver are supplied. A removable plexi-glass plate ;ermits observation during operation.

Chamber Size	Motor Cap.
40 x 25 mm.	0.25 H.P.
65 x 25 mm.	0.25 H.P.
100 x 50 mm.	0.50 H.P.

Lab Willey Grinder Heavy Duty with 1.00 H.P. Motor (Automatic)Modified Model