


BALL MILL

Laboratory Type

www.acmastechnologies.com | www.acmas.in



OBROMAX

ACMAS Technologies Inc.



CE An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603
Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com
www.acmastechnologies.com | www.acmas.in

ACMAS Technologies Inc.

CE An ISO 9001:2000 Company

OBROMAX



BALL MILL Laboratory Type

BALL MILL Laboratory Type

Electrically operated heaving capacity from a few grams, to 1 kg. or 2Kg or 5 Kg. Fitted with a F.H.P. motor cap. up to 2 Kg. jar and with 1/4 H.P. Heavy-duty motor for 5 Kg. jar maximum speed 80 R.P.M. supplied with steel balls for mixing and grinding particularly soft powder. Operates on 220 Volts AC. mains.

- a) 1kg jar
- b) 2kg jar
- c) 5kg Jar.

Capacity

Optional Accessories: Automatic 0-60 minute "German" Timer