

# Constant Bath

Model No: ATI-139



# Constant Bath

Model No: ATI-139



We manufacture and export rectangular, constant temperature water baths also known as double walled construction. Its Inner chamber is made of stainless steel and outer body made of mild steel duly powder coated. Temperature range 5°C above ambient to 90°C.

fitted with circular pump with stirrer.

## Salient features

- CE marked
- Laboratory calibrated
- Robust construction
- Reliable & versatile range
- Energy efficient
- Eco-friendly
- Calibrated equipment
- Controlled environmental conditions
- Completely customized

## Construction & Specification

- Highly recommended where water is required to be circulated to some external source through and F.H.P pump.
- Fitted with stirring with 1/20 H.P. Motor regulator.
- Fitted with immersion type heater.
- Temperature Range is from ambient to 90 deg C with an accuracy of  $\pm 1$  deg C.
- Double walled, Inner made of Stainless Steel, outer made of heavy duty Mild Steel duly powder coated.
- Properly insulated.
- Covered with perforated S.S.tray.
- Temperature is controlled with thermostat.
- Temperature control knob is graduated in degree centigrade.
- Provided with On/Off switch, indicator, cord plug, etc.
- Suitable to work on 220 V single phase 50 Hz AC.

## Technical Specification

Model no.	ATI-139
Temperature range	Ambient to 90 deg C
Temperature accuracy	$\pm 1$ deg C
Controller	(Hydraulic type) Thermostatically Control
Power supply	220 V single phase 50 Hz AC

## Working Chamber

Dimension	
Inside	8" X 9" X 10.5"
Outside(Cover)	10" X 7"

## Optional Accessories:

- Digital temperature indicator cum controller



## ACMAS Technologies Inc.

 An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA

Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603 | Email: acmastechnologies@gmail.com | obromax.delhi@gmail.com

www.acmastechnologies.com | www.acmas.in