

# Rotary Shaker (Horizontal)

# Rotary Shaker (Horizontal)



**OBROMAX**

**ACMAS Technologies Inc.**

CE | An ISO 9001:2000 Company

Sales Office: 312-313, Vardhman Capital Mall, L.S.C. 10, Gulabi Bagh, Delhi-110052, INDIA  
Tel: +91-11-23654603 (M) +91-9313971681, 9311039044, 9555351619 Fax: +91-11-23654603  
Email: [acmastechnologies@gmail.com](mailto:acmastechnologies@gmail.com) | [obromax.delhi@gmail.com](mailto:obromax.delhi@gmail.com)  
[www.acmastechnologies.com](http://www.acmastechnologies.com) | [www.acmas.in](http://www.acmas.in)

**ACMAS Technologies Inc.**

CE | An ISO 9001:2000 Company

**OBROMAX**

Obromax introduced world's best Rotator Shaker is used in blood banks and hematology laboratories for blood typing, cross matching and mixing of blood donor with anti-coagulant, V.D.R.L. Tests and other tests.

**SALIENT FEATURES**

- Operates at 130/180 RPM (as desired) .
- Rotates Specimens in horizontal plane in ¾" circle
- Manufacturer, Exporter and Supplier of Water Bath Shakers
- 30 X 30 cm platform can accommodate slides, blood bottles, flasks & beakers by use of spring bottle holder
- Constructed of all steel finished with enamel paint
- Turned up edges & rubber sheeting on platform to prevent slipping of specimen containers
- Platform rigidly supported by steel legs for smooth & silent rotation and heavy load can be easily be supported.
- Shaker mounted on four rubber suction feet to prevent creeping
- 1/35 H.P motor operates on 220/230 volts AC only Supplied with brass ring maker for 12 rings Without thermometer
- Adjustable 0 to 30 minute timer with arrangement for continuous operation (optional)

**Specifications:**

- Operates at 130/180 RPM (as desired)
- Rotates Specimens in horizontal plane in ¾" circle
- Turned up edges & rubber sheeting on platform to prevent slipping of specimen containers
- 1/35 H.P motor operates on 220/230 volts AC only
- Platform rigidly supported by steel legs for smooth & silent rotation and heavy load can be easily be supported Supplied with brass ring maker for 12 rings
- Adjustable 0 to 30 minute timer with arrangement for continuous operation



**Available in following capacities**

Platform Size	Holding Capacity	HP (Motor)
45 cm × 45 cm	25 flask of 100 mL	0.25
45 cm × 45 cm	36 flask of 50 mL	0.25
45 cm × 45 cm	49 flask of 25 mL	0.25
45 cm × 45 cm	16 flask of 250 mL	0.25
55 cm × 55 cm	25 flask of 250 mL	0.25
65 cm × 65 cm	36 flask of 250 mL	0.50
100 cm × 100 cm	81 flask of 250 mL	1.00
105 cm × 105 cm	81 flask of 500 mL	1.00
105 cm × 105 cm	49 flask of 1000 mL	1.00